

Medi-Stand™ 5 Star



Mobile equipment mounting solution

“Ergonomic design,,

“Flexible,,

“Easy-to-clean,,

“Highly manoeuvrable,,

“Light and Sturdy,,

# Highlights



## Sturdy & Elegant design

Cast aluminium base provides a rigid and steady base for the Medi-Stand™. Its 700mm diameter provides good performance against 'topple' even when loaded with up to 40kg of equipment.

Its elegant form has been designed to be aesthetically pleasing and compliment any medical environment.



## Accessories


Having listened to what medical professionals want, our Medi-Stand™ can be supplied with a standard 25mm x 10mm mounting rail, adjustable drip pole, storage basket and cable tidy hook.



## Adaptable mounting solution

Incorporating the GCX standard mounting section our Medi-Stand™ will accept a wide range of legacy mounting brackets.

Plastic infills keep these mounting channels covered and easy to clean when not in use.



## Handle

Ergonomic aluminium handle is strong and comfortable to use. Using the GCX mounting channel the handle can be positioned and set at a height to suit the operator and any other additional accessories.



## Equipment Mounting plate

The top mounting plate measures 260 x 334mm and is made of sturdy 3mm Mild Steel. Able to support up to 40kg of monitoring equipment, the Medi-Stand™ 5 Star offers great flexibility to the user.

For more information or to discuss your specific requirements please call us on 01737 378 800